

I. Development of a HPLC method for characterization of Kreb's cycle Acids in Tobaccos. (NO derivatization, simple prep.)

Approaches

- Reverse Phase pH=2 10.5% KH_2PO_4 Buffer System.

1. Column effect 2. Linking 2 RCM

- Reverse phase Buffer with organic modifier **optimize condition**

- Ion Exclusion: Interaction Aminex A-9 0.01% H_2SO_4
Aminex HPX-87

- Reverse phase-Ion pair Tetra butyl ammonium phosphate Buffer with ACN

BEST COPY

plus Hexadecyltrimethyl Ammonium chloride

II. Evaluation of New ammoniated /denicotinized Tobacco through Chemical Analysis and Consumer Testing

- G30 FC Tobacco - Ammoniated by GTA process and denicotinized by JDF process.
- Monitor physical, routine and DOF analysis and elucidate the chemical changes caused by processing and storage ageing.
1. Chemical Analysis 2. Amino Sugar, Amadori, DOF

III. Support Technical R&D on Optimization of G7 blend with screening design of ~~important components~~ to find which variables are important. 15 samples