

VAX 26557

Parameter: Total Residual Chlorine
 Method: Amperometric Titration (Direct)

	Yes	No
1. Is sample analyzed immediately after collection?	<input checked="" type="checkbox"/>	<input type="checkbox"/>
2. Is PAO normality 0.00564N?	<input checked="" type="checkbox"/>	<input type="checkbox"/>
3. Are reagents free of contamination or growths?	<input checked="" type="checkbox"/>	<input type="checkbox"/>
4. Is KI solution discarded when it turns yellow?	<input checked="" type="checkbox"/>	<input type="checkbox"/>
5. Is pH 4 buffer solution acetate buffer?	<input checked="" type="checkbox"/>	<input type="checkbox"/>
6. Are reagents within their indicated shelf lives?	<input checked="" type="checkbox"/> 70 shelf life in bottle	<input type="checkbox"/>
7. Is sample volume 200 ml (for chlorine residual up to 2 mg/l) or 100 ml or proportionately less (for chlorine residuals in excess of 2 mg/l)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>
8. Is at least 1 ml KI solution added?	<input checked="" type="checkbox"/>	<input type="checkbox"/>
9. Is at least 1 ml acetate buffer added after KI solution?	<input checked="" type="checkbox"/>	<input type="checkbox"/>
10. Is titrant added in progressively smaller increments until all needle movement ceases?	<input checked="" type="checkbox"/>	<input type="checkbox"/>
11. Is last increment of titrant that causes no needle response subtracted from final volume?	<input checked="" type="checkbox"/>	<input type="checkbox"/>
12. Is the sample value calculated correctly? $[TRC (mg/l) = \frac{A \times 200}{ml \text{ of sample}} \text{ where } A = ml \text{ PAO used}]$	<input checked="" type="checkbox"/>	<input type="checkbox"/>

Comments: SAMPLE IS TAKEN AT THIS TIME (X) ALL REAGENTS IN ACCORDANCE TO THEIR RESPECTS.

Rating: ☒ Satisfactory ☐ Qualified ☐ Unsatisfactory

Problems:

2024407095